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





Rx Only	MD	!	⊗	ENGLISH	Instructions for Use
				Intended Use: This product is intended for use in orthodontic treatment of Class II malocclusions.	
				Product Description: Forsus E22 module is a single-use, non-sterile stainless steel device, designed for use in conjunction with orthodontic brackets, bands, tubes, and archwire products for the treatment of orthodontic malocclusions.	
				Indications for Use: This product is indicated for use in orthodontic treatment of Class II malocclusions.	
				Warning: This product contains nickel and/or chromium. A small percentage of the population is known to be allergic to nickel and/or chromium. If an allergic reaction occurs, direct patient to consult a physician.	
				The Forsus E22 Module requires the use of a buccal tube with an occlusal headgear tube, and is compatible with all such buccal tubes measuring 3.6 mm-4.3 mm in length.	
				Intended Patient Population: Patients requiring orthodontic treatment unless the patient's condition limits the use. The safety and performance of the product has not been clinically evaluated in pregnant or lactating women.	
				Contraindications: None known	
				Intended Users: Educated dental professionals who have theoretical and practical knowledge on usage of dental products.	
				Clinical Benefit: Correct malocclusions when used as part of an orthodontic treatment system.	
				Precaution: Use of recycled product also may compromise material integrity, bond strength, patient safety and overall product performance.	
				Potential Side-Effects: Local irritations, or local or systemic allergic reactions.	
				Residual Risks: Improper placement or breakage can lead to swallowing or aspiration of the product or parts of it.	

Forsus™ Solventum

Instructions For Use	Instrukcja używania
Mode d'emploi	Használati útmutató
Gebrauchsanweisung	Návod k použití
Istruzioni per l'uso	Návod na použitie
Instrucciones de uso	Navodila za uporabo
Gebruiksaanwijzing	Kasutusjuhend
Bruksanvisning	Lietošanas instrukcijas
Brugsanvisning	Naudojimo instrukcija
Bruksanvisning	Instrukcioni de utilizare
Käyttöohjeet	Upute za uporabu
Instruções de uso	Инструкções за употреба
Οδηγίες Χρήσης	Kullanma Talimatları

Solventum

Solventum Orthodontics Corporation
2724 South Peck Road
Monrovia, CA 91016 USA
Solventum.com

Solventum Germany GmbH Edisonstrasse 6 59174 Kamen, Germany Solventum.com	KCI Medical GmbH Eggstrasse 91 8803 Rueschlikon Switzerland		

UK Responsible Person:
KCI Medical Limited, Building 47
Charnwood Campus, 10 Bakewell Road
Loughborough, LE11 5RB, United Kingdom

Issue Date: 2025-08
© Solventum, 2025.
REF 013-634 0272061.1

ENGLISH	Instructions for Use
----------------	-----------------------------

Intended Use: This product is intended for use in orthodontic treatment of Class II malocclusions.

Product Description: Forsus E22 module is a single-use, non-sterile stainless steel device, designed for use in conjunction with orthodontic brackets, bands, tubes, and archwire products for the treatment of orthodontic malocclusions.

Indications for Use: This product is indicated for use in orthodontic treatment of Class II malocclusions.

Warning: This product contains nickel and/or chromium. A small percentage of the population is known to be allergic to nickel and/or chromium. If an allergic reaction occurs, direct patient to consult a physician.

The Forsus E22 Module requires the use of a buccal tube with an occlusal headgear tube, and is compatible with all such buccal tubes measuring 3.6 mm-4.3 mm in length.

Intended Patient Population: Patients requiring orthodontic treatment unless the patient's condition limits the use. The safety and performance of the product has not been clinically evaluated in pregnant or lactating women.

Contraindications: None known

Intended Users: Educated dental professionals who have theoretical and practical knowledge on usage of dental products.

Clinical Benefit: Correct malocclusions when used as part of an orthodontic treatment system.

Precaution: Use of recycled product also may compromise material integrity, bond strength, patient safety and overall product performance.

Potential Side-Effects: Local irritations, or local or systemic allergic reactions.

Residual Risks: Improper placement or breakage can lead to swallowing or aspiration of the product or parts of it.

The following risks may occur: tissue abrasion or laceration, tooth sensitivity, tooth damage or fracture, pain or discomfort, occlusal changes or temporal mandibular joint discomfort.

Handling Before and/or After Use: This product is for single use only. Properly dispose of used or contaminated product to avoid health risks of improper handling.

SDS and Disposal: See the Solventum Safety Data Sheet (SDS)/Safety Information Sheet (SIS), available at www.Solventum.com or contact your local subsidiary, for disposal and ingredient information.

Dispose of contents or container in accordance with the applicable regulations.

Please pay special attention to the disposal of contaminated waste, waste with sharps to avoid health risks because of improper handling.

Customer Information: Please report a serious incident that has occurred in relation to the device to Solventum and the local competent authority (EU) or local regulatory authority.

Symbols Glossary: For symbols glossary, see eIFU.Solventum.com

Forsus E22 Module Installation

Holding the mesial portion of the Forsus™ Fatigue Resistant Device E22 Module with a pair of Weingart Utility Pliers (Figure 1), insert the Forsus E22 Module into the headgear tube from mesial to distal until it clicks into place. (Figure 2) **(Note: Right and Left configuration denoted by "R" or "L" on each module).**

Cinch archwire to secure lower arch teeth positioning. (Figure 3 and 4) For increased anchorage on lower arch, use second molar attachment if appropriate. **(Note: Use full size rectangular stainless steel wire on both arches to further secure teeth positioning).**

In order to avoid elastic ligature breakage and to secure cuspid bracket to archwire, use stainless steel ligature tie. (Figure 5)

Push Rod Installation

Measure from distal end of upper molar tube to the chosen mandibular stop while in centric occlusion. Use Measurement Guide (REF 807-014) as shown in picture and note appropriate push rod size (in this example, note 32 mm). (Figure 6)

Select a push rod from the available sizes to fully compress the spring without immediately reposturing the mandible. (Figure 7) **(Note: Right and Left configuration).**

Place push rod loop on archwire (Figure 8), have patient open mouth, compress spring, and insert push rod. If push rod protrudes distal of spring module in centric occlusion, it is too long. (Figure 9)

Crimp push rod mesial end by closing loop around archwire. (Figure 10)

Re-Activation

Reactivation may only be achieved on push rod sizes that have stops: 25 mm, 29 mm, 32 mm, and 35 mm. To re-activate, crimp split ring bushings on push rod distal of stop by compressing the spring as needed (usually 2 or 3 mm at a time to achieve midline correction and advancement). (Figure 11)

Forsus E22 Module Removal

First, with the patient's mouth open wide, compress the spring and remove the spring assembly from the push rod. Then, holding the distal portion of the Forsus E22 Module with a pair of Weingart Utility Pliers (Figure 12), pull the Forsus E22 Module in an occlusal/mesial direction to remove it from the headgear tube.

Push Rod Removal

Grasp the distal end of the push rod with your fingers. Holding the mesial loop in place with your Weingart plier (Figure 13), bend the push rod to open the mesial loop. (Figure 14)